

# Surgical Suture Materials

CATGUT GMBH  
MARKNEUKIRCHEN / Germany

CATGUT GMBH MARKNEUKIRCHEN / Germany		Structure of Thread	Gauge	Colour	Tensile strength	Time of degradation
Absorbable Sutures:					appr.	appr.
<b>MARLIN®</b>	(polyglycolic acid)	monofilam.:10-0 to 9-0 braided/ coated	USP 10-0 to 9-0 USP 8-0 to 6	violet/ undyed	50% 14-16 days 0% 28-35 days	90 days
<b>MARLIN® rapid</b>	(polyglycolic acid)	braided/ coated	USP 8-0 to 2	undyed violet:Gauges 8-0, 7-0	50% 7 days 0% 21 days	42 days
<b>PCL</b>	(PolyCaprolacton-co-Lactid)	monofilament/ coated	USP 8-0 to 2	violet	50% 42-50 days 0% 90-100 days	180-210 days
<b>MARISORB</b>	(Poly-p-dioxanon)	monofilament	USP 7-0 to 2	violet/ undyed	50% 42-50 days 0% 90-100 days	180-210 days
<b>MARICRYL®</b>	Glykolide-ε-Caprolactone- Copolymer	monofilament	USP 6-0 to 1	undyed/ violet	50% 7 days 0% 28 days	90-120 days
Non-absorbable Sutures:						
<b>Marilon®</b>	(Polyvinylidene fluoride)	monofilament	USP 10-0 to 2	blue	-	-
<b>Mariderm</b>	(Copolyamide 6-6,6)	monofilament	USP 6-0 to 2	blue	-	-
<b>Mariderm</b>	(Polyamide 6,6)	monofilament	USP 11-0 to 6-0	black	-	-
<b>Maprolen®</b>	(Polypropylene)	monofilament	USP 10-0 to 1	blue	-	-
<b>Polyester-S green</b>	(Polyethylene terephthalate)	braided/ coated	USP 6-0 to 5	green	-	-
<b>Polyester white</b>	(Polyethylene terephthalate)	braided	USP 7-0 to 7	white	-	-
<b>Greenfil</b>	Polyester (Polyethylene terephthalate)	braided	USP 6-0 to 4	green	-	-
<b>Silk-S</b>		braided/ coated	USP 8-0 to 6	black	-	-
<b>Steel wire</b>		monofilament multifilament	USP 4-0 to 2 USP 3-0 to 2	natural natural	- -	- -